

BINDER BENEFITS:

Improves tablet quality: hardness and weight control Improves particle size control and reduces fines Improves powder compaction in capsules Reduces amount of other excipients required Reduces defects: friability, breakage and capping Reduces compression force required

APPLICATIONS:

Tablets and capsules

ALTERNATIVE TO:

Polyvinylpyrrolidone (PVP) Hydroxypropyl cellulose (HPC) Pregelatinized Starch PolyethyleneGlycol (PEG) Hydroxypropyl Methylcellulose (HPMC)

Gelatin

CLEAN LABEL RECOMMENDATION:

Gum Fiber Blend



Nu-BIND® USAGE RATES

Recommended to start at 1:1 ratio

Binder - 2-5%

INSTRUCTIONS:

For use as a dry binder, add Nu-BIND® to the blend during the final stage of mixing. Blend until homogeneous (approximately 10 minutes). If additional binding capability is needed, hydrate a portion of the Nu-BIND® for use as a liquid addition.

May be used as a wet granulating agent as well.

INGREDIENTS:

Guar Gum, Gum Arabic, Agave Fiber, Rice Hulls and Agave Syrup

(Photo of Nu-BIND®)







